

Servicing the Philips BV-Pulsera C-Arm



RADIOLOGICAL SERVICE TRAINING INSTITUTE

Introduction

This course is designed to provide the advanced service professional with the skills and knowledge to maintain the Philips BV-Pulsera at the highest state of readiness. All adjustments will be discussed to establish optimum performance criteria. Theory and hands-on sessions will develop the skills necessary to troubleshoot system failures and restore it to operation.

Prerequisites

To attend this course, the service professional must have good fundamental knowledge and understanding of the principles gained through attendance at our Phase I, Phase II, and Phase III X-ray courses or equivalent field experience.

Objectives

At the completion of this course participants will be able to:

- Operate the Philips BV-Pulsera
- Identify all systems, subsystems and components of the Philips BV-Pulsera Verify power supplies for accuracy and function
- Service and calibrate system batteries and charger circuits
- Utilize all communication interfaces to calibrate and evaluate the systems
- Evaluate the performance of the X-ray generator, imaging and workstation sections of each system

- Calibrate and adjust all components of the X-ray generator, imaging chain and workstation
- Utilize all diagnostic indicators to troubleshoot system failures
- Restore the system to proper functional state following a system failure
- Evaluate and repair mechanical systems
- Load system software

Course Outline

Day 1

- Introduction
 - o Course objectives
 - o System
 - Major components
 - Configurations
 - Documentation
 - o System operation
 - o C-arm controls
 - X-ray subsystem
 - I.I., collimator and CCD camera controls
 - Mechanical systems
 - o Workstation controls
- Physical layout and component identification
- AC power distribution
- DC power distribution
- Batteries and charger
- Lab Activities
 - o System operation
 - Fluoroscopic modes
 - Low dose

- High dose
- Boost
- Radiographic mode
- Patient data input
- Recall stored images
- Collimator controls
- TV/ II controls
- o Physical layout and component ID
 - Covers and panels
 - Power supplies
 - Circuit boards
 - Battery removal
 - X-Ray tube removal and installation
 - I.I. removal and replacement
 - Mechanical systems
- o Power supply verifications
- o Battery charger calibration
- o Battery charger test points and waveforms

Day 2

- System communications
- Interlocks
- Calibration software interface
 - o Direct connect procedure
 - o Using BV-SCOPE
- X-ray generator
 - o Stator power and control
 - o Pre-charge
 - o X-ray on, x-ray disable
 - o High voltage control
 - o Filament/ mA control
- Lab Activities

NOTE:

Due to copyright laws, students are required to purchase and bring to class a copy of OEC's Service Documentation CD.

9900 P/N 5325033-2EN, Rev001

To order this CD, Please call OEC/GE (800)874-7378.

Servicing the Philips BV-Pulsera C-Arm



RADIOLOGICAL SERVICE TRAINING INSTITUTE

- o Verify ARCNET communication
- o Initialize calibration modes for 9900
- o RUS communication
 - Calibration screens
 - Logging
 - Calibration process
- o Enable system monitors
- o Connect external keyboard
- o Initiate service mode for BV-Pulsera
- o Verify stator operation
- o Pre-charge test
- o Verify x-ray enable signals
- o Calibrate x-ray generator
- o High voltage test points and waveforms
- o Filament/mA control test points and waveforms
- o Max "R" adjustment

Day 3

- Imaging system components
 - o X-ray tube
 - Central ray adjustment
 - Filtration
 - o Image/Fluoro functions control PCB's
 - Collimator control
 - Collimator iris size and center
 - Semi-transparent leafs
 - Width
 - Rotation
 - Image intensifier
 - Size control
 - Focusing
 - Camera
 - Centering

- Focus
- TV camera iris
 - Thermoelectric cooling
- Dose/ brightness control
 - o Video path
 - o Brightness control processing
 - o Iris adjustment
- Lab Activities
 - o Central ray adjustment
 - o Collimator centering
 - o Collimator size tracking calibration
 - o Verify and adjust II viewed field
 - o I.I. focus adjustment
 - o TV camera focus
 - o TV camera center adjustment
 - o TV camera rotation adjustment
 - o TEC verification and adjustment
 - o TV camera iris adjustment

Day 4

- Image display
 - o Image processor
 - Image manipulation
 - Window/level
 - Subtraction
 - Noise suppression
 - Motion artifact suppression
 - o Video distribution board
 - Video input
 - High resolution video output
 - Standard resolution video output
 - o Monitors
 - o Touch screen
- Image storage
 - o Single disk
- Mechanical systems
 - o Flip flop

- o Orbital motion
- o Wig-wag motion
- o Horizontal cross arm motion
- o L-Arm
- o Vertical lift
- o Steering and breaking
- Lab Activities
 - o Monitor adjustments
 - o Image centering
 - o Vertical lift drive tests
 - o Wig-wag adjustment
 - o Mechanical evaluation

Day 5

- Diagnostics
 - o Error messages
 - o LED displays
 - o Bar graphs
 - o Seven segment displays
 - o Monitor "debug" screens
 - o Status monitor
 - o Status/error logs
- System software installation
- Lab Activities
 - o Evaluate for diagnosis
 - LED functions
 - Bar graphs
 - Seven segment displays
 - o Use service mode to verify system operation
 - o View status logs
 - o View error logs
 - o Reload system software
 - o System troubleshooting
 - o System review
 - o Course evaluation